
**sylc™ Propy Therapy**  
Therapeutic Propy Powder

POWERED BY **NOVAMIN**

## DIRECTIONS FOR USE

### Calcium Sodium Phosphosilicate cleaning and desensitization powder

**CAUTION:** Federal law restricts this device to be used by or on the order of a licensed dentist.

**DESCRIPTION:** Sylc is designed to provide a prophylaxis treatment for dental enamel and to physically occlude dentin tubules for the management of sensitive teeth. Sylc is a dry powder (calcium sodium phosphosilicate), composed of elements that occur naturally in the body's hard tissues (Ca, Na, Si, P, and O). When exposed to an aqueous environment, Sylc undergoes a rapid surface reaction, allowing it to physically adhere to exposed dentin and to physically occlude tubules. Within a short period of time, essentially all of the Sylc reacts to form hydroxycarbonate apatite (HCA), which is chemically similar to natural tooth mineral. Sylc removes plaque, soft deposits, and surface stains. It is also recommended for the following procedures:

- Removal of extrinsic stains, e.g. coffee, tea and tobacco
- Treatment of hypersensitivity
- Surface preparation prior to bonding and sealants
- Before and after scaling and root planning procedure to reduce sensitivity
- Prophylaxis of orthodontic patients
- Plaque removal prior to fluoride treatments
- Pre and post whitening treatments

**INDICATIONS:** Sylc is intended for cleaning and polishing procedures as part of a dental prophylaxis treatment. The OSspray cleaning powder is also intended to provide relief from tooth sensitivity during the prophylaxis treatment. Sylc is intended to be used in OSspray recommended air polishing equipment or handheld devices.

**INGREDIENTS:** Sylc particulate (calcium sodium phosphosilicate).

**CONTRAINDICATIONS:** Sylc should not be considered for patients with a known allergy to silica.

### PRECAUTIONS:

- During the preparation process and application, it is recommended that both clinician and patient wear protective eyewear.
- If any material gets into the eyes, immediately flush with water.
- Delivering Sylc with air polishing equipment at high air pressures (psi), may cause injury to the gum tissues or an emphysema caused by the introduction of air into the soft tissue spaces.  
It is strongly recommended not to direct Sylc at the gum tissue or into the gingival sulcus.
- Avoid directing Sylc powder toward restorations, crowns and bridgework as this could damage the restoration.
- Avoid directing Sylc powder on orthodontic brackets as this could dull the metal.
- Empty the powder chamber daily, especially in areas of the country with high humidity to avoid clogging.
- Propy tip lifetime will be optimized when working no higher than the recommended operating air pressure.
- Safety and effectiveness in pregnant women and children have not been established.

### USAGE INSTRUCTIONS:

1. Sylc is intended to be used in OSspray recommended air polishing equipment or handheld devices.
2. The Sylc bottle is opened by turning the lid counterclockwise.
3. Pour Sylc powder into the delivery system (chamber) at the maximum level recommended by equipment manufacturer.
4. OSspray recommends Sylc be applied at an air pressure range of 40-46 psi or approximately 2.8-3.2 bar. Equipment should be set to a minimum powder flow and powder uptake.
5. If increased cleaning effect is required, increase flow and powder settings accordingly.
6. Following treatment, the patient should not eat or drink for at least 60 minutes.
7. Sylc can be used on patients who are on a sodium-restricted diet.

### STABILITY:

- Do not use if bottle or sachet is opened or damaged.
- Do not use after expiration date.
- Store in a DRY place at room temperature, between 15-25 °C and 60-75 °F.

### HOW SUPPLIED:

Sylc is supplied in bottles or unit dose packs.

### TECHNICAL CHANGES:

OSspray reserves the right to make changes to Sylc usage instructions when required in the course of ongoing technical and scientific developments.

### MANUFACTURER:

OSspray Ltd, Hodgkin Building, King's College London, Guy's Hospital Campus, London, SE1 1UL, United Kingdom

OSspray Inc, 301 Pleasant Street, Abbottstown, PA 17301, U.S.A. 866.238.9902

Sylc® is a registered trademark of OSspray, Ltd. Novamin® is a registered trademark of GlaxoSmithKline plc.

**OSSPRAY**

Building Stronger Smiles. Naturally.

[www.osspray.com](http://www.osspray.com)

© 2011, OSspray, Ltd. All rights reserved.  
PM-0024 Rev.0004 04/2011 ENGLISH

CE

**sylc™ Propy  
Therapy**  
Therapeutic Propy Powder

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© 2011, OSsray, Ltd. All rights reserved.